



Datasheet

HCA-S-400M-SI

400 MHz Photoreceiver with Si PIN Photodiode



The picture shows the HCA-S-400M-SI-FS with free space input. The photoreceiver will be delivered without post holder and post.

Features	<ul style="list-style-type: none"> • Si PIN Detector, 0.8 mm Active Diameter • Spectral Range 320 ... 1000 nm • Bandwidth DC ... 400 MHz • Amplifier Transimpedance (Gain) 5.0×10^3 V/A • Max. Conversion Gain 2.7×10^3 V/W @ 800 nm 																															
Applications	<ul style="list-style-type: none"> • Spectroscopy • Fast Pulse and Transient Measurements • Optical Triggering • Optical Front-End for Oscilloscopes and A/D Converters 																															
Specifications	<table border="0"> <tr> <td></td> <td><i>Test Conditions</i></td> <td><i>$V_s = \pm 15$ V, $T_a = 25^\circ$C</i></td> </tr> <tr> <td rowspan="2">Gain</td> <td>Transimpedance</td> <td>5.0×10^3 V/A (@ 50 Ω load)</td> </tr> <tr> <td>Max. Conversion Gain</td> <td>2.7×10^3 V/W (@ 800 nm)</td> </tr> <tr> <td rowspan="4">Frequency Response</td> <td>Lower Cut-Off Frequency</td> <td>DC</td> </tr> <tr> <td>Upper Cut-Off Frequency (-3 dB)</td> <td>400 MHz (± 10 %)</td> </tr> <tr> <td>Rise/Fall Time (10% - 90%)</td> <td>1.0 ns</td> </tr> <tr> <td>Gain Flatness</td> <td>± 1 dB</td> </tr> <tr> <td rowspan="3">Detector</td> <td>Detector Material</td> <td>Si PIN photodiode</td> </tr> <tr> <td>Active Area</td> <td>\varnothing 0.8 mm</td> </tr> <tr> <td>Spectral Response</td> <td>320 ... 1000 nm</td> </tr> <tr> <td rowspan="3">Input</td> <td>Input Offset Compensation Range</td> <td>± 200 μA adjustable by offset trimpot</td> </tr> <tr> <td>Optical Saturation Power</td> <td>400 μW (for linear amplification, @ 800 nm)</td> </tr> <tr> <td>Min. NEP</td> <td>40 pW/\sqrtHz (@ 800 nm, 100 MHz)</td> </tr> </table>		<i>Test Conditions</i>	<i>$V_s = \pm 15$ V, $T_a = 25^\circ$C</i>	Gain	Transimpedance	5.0×10^3 V/A (@ 50 Ω load)	Max. Conversion Gain	2.7×10^3 V/W (@ 800 nm)	Frequency Response	Lower Cut-Off Frequency	DC	Upper Cut-Off Frequency (-3 dB)	400 MHz (± 10 %)	Rise/Fall Time (10% - 90%)	1.0 ns	Gain Flatness	± 1 dB	Detector	Detector Material	Si PIN photodiode	Active Area	\varnothing 0.8 mm	Spectral Response	320 ... 1000 nm	Input	Input Offset Compensation Range	± 200 μ A adjustable by offset trimpot	Optical Saturation Power	400 μ W (for linear amplification, @ 800 nm)	Min. NEP	40 pW/ \sqrt Hz (@ 800 nm, 100 MHz)
	<i>Test Conditions</i>	<i>$V_s = \pm 15$ V, $T_a = 25^\circ$C</i>																														
Gain	Transimpedance	5.0×10^3 V/A (@ 50 Ω load)																														
	Max. Conversion Gain	2.7×10^3 V/W (@ 800 nm)																														
Frequency Response	Lower Cut-Off Frequency	DC																														
	Upper Cut-Off Frequency (-3 dB)	400 MHz (± 10 %)																														
	Rise/Fall Time (10% - 90%)	1.0 ns																														
	Gain Flatness	± 1 dB																														
Detector	Detector Material	Si PIN photodiode																														
	Active Area	\varnothing 0.8 mm																														
	Spectral Response	320 ... 1000 nm																														
Input	Input Offset Compensation Range	± 200 μ A adjustable by offset trimpot																														
	Optical Saturation Power	400 μ W (for linear amplification, @ 800 nm)																														
	Min. NEP	40 pW/ \sqrt Hz (@ 800 nm, 100 MHz)																														


400 MHz Photoreceiver with Si PIN Photodiode

Specifications (continued)																										
Output	Output Voltage Range	± 1.0 V (@ 50 Ω load) for linear operation and low harmonic distortion																								
	Max. Output Voltage Range	± 1.5 V (@ 50 Ω load)																								
	Output Impedance	50 Ω (designed for 50 Ω load)																								
	Output Noise	ca. 20 mV peak-peak or 3 mV rms (@ 50 Ω load, no signal on detector)																								
Power Supply	Supply Voltage	± 15 V																								
	Supply Current	± 55 mA typ. (depends on operating conditions, recommended power supply capability minimum ± 150 mA)																								
Case	Weight	210 g (0.5 lbs)																								
	Material	AlMg4.5Mn, nickel-plated																								
Temperature Range	Storage Temperature	- 40 ... + 100 °C																								
	Operating Temperature	0 ... + 60 °C																								
Absolute Maximum Ratings	Optical Input Power	20 mW																								
	Power Supply Voltage	± 22 V																								
Spectral Response	<table border="1"> <caption>Spectral Response Data (Estimated)</caption> <thead> <tr> <th>Wavelength [nm]</th> <th>Photo Sensitivity [A/W]</th> </tr> </thead> <tbody> <tr><td>300</td><td>0.10</td></tr> <tr><td>350</td><td>0.15</td></tr> <tr><td>400</td><td>0.12</td></tr> <tr><td>500</td><td>0.20</td></tr> <tr><td>600</td><td>0.35</td></tr> <tr><td>700</td><td>0.48</td></tr> <tr><td>800</td><td>0.55</td></tr> <tr><td>850</td><td>0.55</td></tr> <tr><td>900</td><td>0.50</td></tr> <tr><td>1000</td><td>0.20</td></tr> <tr><td>1050</td><td>0.10</td></tr> </tbody> </table>		Wavelength [nm]	Photo Sensitivity [A/W]	300	0.10	350	0.15	400	0.12	500	0.20	600	0.35	700	0.48	800	0.55	850	0.55	900	0.50	1000	0.20	1050	0.10
Wavelength [nm]	Photo Sensitivity [A/W]																									
300	0.10																									
350	0.15																									
400	0.12																									
500	0.20																									
600	0.35																									
700	0.48																									
800	0.55																									
850	0.55																									
900	0.50																									
1000	0.20																									
1050	0.10																									
Connectors	Input	HCA-S-400M-SI-FS 25 mm round flange for free space applications HCA-S-400M-SI-FC FC fiber optic receptacle HCA-S-400M-SI-SMA SMA fiber optic receptacle																								
	Output	BNC																								
	Power Supply	LEMO series 1S, 3-pin fixed socket Pin 1: + 15V Pin 2: - 15V Pin 3: GND																								


400 MHz Photoreceiver with Si PIN Photodiode

Typical Performance Characteristics

Frequency Response



Noise Spectrum



Note: Spectral noise data is measured at the amplifier output with no signal on the photodiode. To determine the spectral input noise divide the measured output noise by the amplifier conversion gain.


Conversion gain (V/W) = amplifier gain (5,000 V/A) x photo sensitivity (A/W).

Marker	Frequency	Output Noise	Resulting Input Noise (NEP)
1	100 MHz	109 nV/√Hz	40 pW/√Hz (@ 800 nm)
2	380 MHz	183 nV/√Hz	68 pW/√Hz (@ 800 nm)


400 MHz Photoreceiver with Si PIN Photodiode

Typical Performance Characteristics (continued)

Pulse Response to Square Wave Input Signal (with 16 times averaging)




Large Signal Response output signal for 100 MHz, 370 μ W modulated optical input signal (with 4 times averaging)




400 MHz Photoreceiver with Si PIN Photodiode

Typical Performance Characteristics (continued)

Small Signal Response
output signal for 3.7 μ W modulated optical input signal, 1 MHz square wave (without (top) and with 64 times averaging (bottom))





CH1 INPUT
DC50 Full
10.0mV/div
1:1




CH1 INPUT
DC50 Full
10.0mV/div
1:1

400 MHz Photoreceiver with Si PIN Photodiode

Available Models	HCA-S-400M-SI-FS HCA-S-400M-SI-FC HCA-S-400M-SI-SMA HCA-S	free space input FC fiber optic receptacle SMA fiber optic receptacle customized versions available on request
------------------	--	---

Dimensions	<p>HCA-S-400M-SI-FS</p>  <p>all measures in mm unless otherwise noted DZ-HCA-S-FS_R2</p> <p>HCA-S-400M-SI-FC</p>  <p>all measures in mm unless otherwise noted DZ-HCA-S-FC_R4</p>	
------------	--	--

400 MHz Photoreceiver with Si PIN Photodiode



FEMTO Messtechnik GmbH
 Klosterstr. 64
 D-10179 Berlin · Germany
 Tel.: +49-(0)30-280 4711-0
 Fax: +49-(0)30-280 4711-11
 e-mail: info@femto.de
 http://www.femto.de

Specifications are subject to change without notice. Information furnished herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights granted by implication or otherwise under any patent rights of FEMTO Messtechnik GmbH. Product names mentioned may also be trademarks used here for identification purposes only.
 © by FEMTO Messtechnik GmbH
 Printed in Germany